



## Prymax Angi Inflation Device


### DEVICE DESCRIPTION:

An inflation device specifically designed for dilation balloons. This inflation device is a 20 ml disposable device capable of producing a maximum pressure of 30 ATM.

### AVAILABLE VARIANTS:

Straight type, Gun type

### Drawing of product:



### MATERIAL USED:



Poly Carbonate, Acrylonitrile-butadiene-styrene, Polyurethane

### INTENDED PURPOSE:

INTENDED USE: Inflation Device is used to inflate & deflate the balloon and measure the pressure within the balloon during the PTCA procedure

INTENDED USER: The device is intended to be used only by a qualified medical professional.

INTENDED PATIENT POPULATION: All patient population irrespective of age and gender.

Prepared by:	Approved By:
	

**INDICATIONS:**

Percutaneous transluminal coronary angioplasty

**CONTRAINDICATIONS:**

- This device is not designed, sold, or intended for use except as indicated.

**CLINICAL BENEFITS & PERFORMANCE CHARACTERISTICS****Clinical Benefits:**

- Inflation devices help minimize risks and ensure patient comfort by providing consistent pressure regulation and accurate measurement.
- Inflation devices help improve procedural efficiency by facilitating proper vessel dilation, accurate stent deployment, or controlled fluid administration.

**PERFORMANCE CHARACTERISTICS:**



Accurate pressure control

**INSTRUCTION FOR USE:**

- Open the product in sterile environment.
- Visually inspect the product & packing for any visual defects. Improper transportation & handling can cause structural/functional damage to the device or packaging.
- The device is designed for maximum application pressure of 30 bars.
- Hold the device and push the buttons forward to deactivate the lock.
- By keeping the lock in deactivated position, pull the plunger handle back to fill the radiography solution into the device.
- Release the button to activate the lock.
- Plunger handle can not be pushed or pulled unless the lock is deactivated
- Remove air by pushing plunger handle.
- Connect the Male Luer of Inflation Device to the balloon catheter. Secure the connection. Follows the indications and instructions for balloon catheter procedure as prescribed by the manufacturer.
- Once the balloon catheter is positioned, screw the plunger handle in clockwise direction to increase the pressure and thus, for inflation of the balloon.
- Screwing the plunger handle in anti-clockwise direction will reduce the pressure and thus, for deflation of the balloon.

**WARNING & PRECAUTIONS:**

- The use of this product is restricted to a qualified doctor or a surgeon trained in the technique.
- Read all the information mentioned in instruction for use. Failure to do so may result in severe patient injury or death. The product should be used according to the instructions for use.
- The product is guaranteed, sterile till the package has not been opened or damaged within the expiry date.

Prepared by:	Approved By:
	

- Do not use if the sterile pack is opened or damaged.
- Sterile, Single use: Do not reuse, reprocess or re-sterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.
- Always follow the manufacturer's directions accompanying the balloon dilation catheter for instructions for use, maximum balloon inflation pressure, precautions, and warnings for that device. Troubleshooting: User Paramedic staff should be aware & able to troubleshoot the certain complications that may occur during the application of the product such as; He/she should be careful while performing the balloon angioplasty to avoid rupturing or dissecting the vessel wall.

**ADVERSE EFFECTS:**

Vascular damage, Air embolism

**RE-USE HAZARDS:**

Transmission of infection from one patient to another.

**PACKAGING:**

The device is packed in Header Pouch, which is further packed in inner carton & master carton.

**STERILIZATION:**

Device is EO sterilized, and sterilization is done in-house.

**STORAGE CONDITION:**

Temperature Limit: 15°C to 40°C.

**DEVICE SHELF LIFE:** 3 Year

**COUNTRY CODE:** 356

**DISPOSAL**

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state, and federal law and regulation.



**REPORTING OF ADVERSE EVENTS:**

In case of any serious adverse event please report the same to [info@prymaxhealthcare.com](mailto:info@prymaxhealthcare.com) & [gloizou@compliancems.com.au](mailto:gloizou@compliancems.com.au).

Note: For product used in EU also report to the competent authority of the Member State in which the user and/or patient is established.

**GOODS RETURN POLICY**

Refer to the company's return good policy.

Prepared by:	Approved By:
	

Please contact the branch office or customer service at [info@prymaxhealthcare.com](mailto:info@prymaxhealthcare.com) or call +91 129 4011818.

#### SYMBOLS USED ON PRODUCT LABEL



Lot number



Do not re-use



Sterilized by Ethylene Oxide



Manufacturing date



Use-by date



Caution



Non-pyrogenic



Do not re-sterilize



Catalogue number



Manufacturer



Storage Temperature Range



Conformité Européenne



Do not use if package is damaged



Consult instructions for use




Single Sterile Barrier System



Medical Device



UDI Carrier



Country of Manufacturer





**Manufactured By:**

**Prymax Healthcare LLP.,**  
53/17, Industrial Area NIT, Faridabad,  
Haryana, INDIA -121001  
[www.Prymaxhealthcare.com](http://www.Prymaxhealthcare.com)



**Compliance Management Solutions Europe LTD**  
2 Bulgaria Str. ,2850 Petrich, Bulgaria  
E-mail: [gloizou@compliancems.com.au](mailto:gloizou@compliancems.com.au)  
Contact No.: **+61 4 33 124266**

Prepared by:	Approved By:
	

Confidential

Proprietary Information of Prymax Health Care LLP  
MN/PAID/001, Rev:02, Date:20/11/2023